

Disinfectant

MILzyme

Triple
Enzyme

Milzyme contains surfactants and triple enzymes, providing a solution for the effective and practical washing of surgical instruments, medical devices and fabrics

HOW TO USE :

To clean the instrument or cloth, first dissolve 35 mL of Milzyme in 3,5 litres of water. Then place the item to be cleaned into the solution and soak for 5-30 minutes. Rinse with clean water several times if necessary, then rub to remove any remaining residue.

ACTIVE INGREDIENTS

Enzyme.....7%
Surfactant.....20%

STORAGE

Store at normal room temperature and dry
Keep out of reach of children
Protect from light

WARNING

In case of contact with eyes and other mucous tissues, rinse immediately with water. External use only

- Does not damage instruments made of metal, plastic, rubber and the like.
- Suitable for cleaning a variety of surgical instruments and medical devices.

5L



Expired :

Production Code :

KEMENKES RI AKD 20903020017

Manufactured by :
CV. Samudra Medika Laboratories
Bandung - Indonesia